

FEATURES


- 10-Channel & 16bands WDRC
- Intelligent noise reduction
- Fuzzy compression detectors
- Static feedback cancellation
- Multi-Memory Selection
- Low level expansion
- Personal fitting kits
- Interactive fitting experience




Technical Data (ANSI S3.22-2003)

Frequency Range	100-6500 Hz	
Peak Gain	60	
Peak Output	125 dB	
HFA Gain	55 dB	
Equivalent Input Noise	<25 dB	
THD at RTG	500 Hz	2%
	800 Hz	1%
	1600 Hz	1%
time constan	Attack	5 ms
	Release	200 ms
Battery type	13	
Battery life	about 210h	

Suggest fitting range



Maximal Output



Maximal Gain

