




Model	Intelligia T ITC			Style	ITC
Fitting Range	Mild to severe hearing loss				
Channels	10				
Bands	15				
Matrix	118/50/X				
Battery	312	Current Consumption	< 1 mA		

FUNCTION	• WDRC • Feedback suppression • Dual time constants • Low Level Expansion • Multi-memory • Noise attenuation algorithm • Adjustable noise reduction • Multi TK fitting	Manual Volume Control	Yes
		Programmable Parameters	• Gain at each channel • TK at each channel • MPO at each channel • MPO TK at each channel • Frequency Shaping • LLE at each channel • NAA 1/2/OFF


ANSI S3.22-1996 Technical Data

Frequency Range		200-7500 Hz
Peak Gain		50 dB
Peak Output		118 dB
Reference Test Gain		40 dB
HFA Output		115 dB
HFA Gain		48 dB
Equivalent Input Noise		<25 dB
THD	Typical	
	500 Hz	1 % < 5%
	800 Hz	1 % < 5%
	1600 Hz	1 % < 5%
SLOW	Attack Time	50 ms
	Release Time	500 ms
FAST	Attack Time	5 ms
	Release Time	20 ms

Fitting Range



OSPL90 Curve



Gain Curve

